

SAFETY DATA SHEET

according to Regulation (EC) No. 1907/2006

Aspirmatic **No Change Service!**

Version 01.07

Revision Date 23.03.2009

Print Date 02.09.2009

Skin contact : Wash off immediately with plenty of water.
Ingestion : Obtain medical attention.

5. FIRE-FIGHTING MEASURES

Special protective equipment for fire-fighters : No special precautions required.
Suitable extinguishing media : Water, Dry powder, Foam, Carbon dioxide (CO₂)
Extinguishing media which shall not be used for safety reasons : No information available.
Specific risk from the substance or the product itself, its combustion products or evolved gases : No special risks to be expected.

6. ACCIDENTAL RELEASE MEASURES

Personal precautions : Increased risk of slipping in the presence of leaked / spilled product.
Environmental precautions : Do not flush into surface water.
Methods for cleaning up : Wipe up with absorbent material (e.g. cloth, fleece). Soak up with inert absorbent material (e.g. sand, silica gel, acid binder, universal binder, sawdust).

7. HANDLING AND STORAGE

Handling

Safe handling advice : Prepare the working solution as given on the label(s) and/or the user instructions.
Advice on protection against fire and explosion : No special protective measures against fire required.

Storage

Requirements for storage areas and containers : Store at room temperature in the original container.
Further information : Keep away from heat. Keep container tightly closed.
Advice on common storage : Keep away from food and drink.

8. EXPOSURE CONTROLS/PERSONAL PROTECTION

Components with workplace control parameters

| Components | CAS-No. | Value | Basis |
|------------|---------|--|----------|
| Ethanol | 64-17-5 | AGW: 960 mg/m ³ AGW: 500 ppm | TRGS 900 |

Personal protective equipment

Hand protection : Splash protection: disposable nitrile rubber gloves e.g. Dermatril (Layer thickness: 0,11 mm) made by KCL or gloves from other manufacturers offering the same protection.

SAFETY DATA SHEET

according to Regulation (EC) No. 1907/2006

Aspirmatic **No Change Service!**

Version 01.07

Revision Date 23.03.2009

Print Date 02.09.2009

Prolonged contact: Nitrile rubber gloves e.g. Camatril (>480 Min., Layer thickness: 0,40 mm) or butyl rubber gloves e.g. Butoject (>480 Min., Layer thickness: 0,70 mm) made by KCL or gloves from other manufacturers offering the same protection.

Eye protection : Safety glasses
Hygiene measures : Keep away from food and drink.
Protective measures : Avoid contact with skin and eyes.

9. PHYSICAL AND CHEMICAL PROPERTIES

Appearance

Form : liquid
Colour : blue
Odour : characteristic

Other data

solidification range : ca. 0 °C
Initial boiling point : ca. 100 °C
Flash point : > 100 °C
Method: ISO 2719
Density : ca.0,99 g/cm³ at 20 °C
Water solubility : at 20 °C
in all proportions
pH : ca.7 at 1.000 g/l (20 °C)

10. STABILITY AND REACTIVITY

Hazardous reactions : None reasonably foreseeable.
Hazardous decomposition products : None reasonably foreseeable.

11. TOXICOLOGICAL INFORMATION

Further information : No data is available on the product itself. The classification was made according to the calculation procedure of the Preparations Directive.

12. ECOLOGICAL INFORMATION

Biodegradability : Readily biodegradable.
Method: OECD 301D / EEC 84/449 C6
Toxicity to bacteria : EC50 520 mg/l

SAFETY DATA SHEET

according to Regulation (EC) No. 1907/2006

Aspirmatic **No Change Service!**

Version 01.07

Revision Date 23.03.2009

Print Date 02.09.2009

Method: OECD 209

Chemical Oxygen Demand (COD) : ca. 2630 mg/l
Test substance: 1% solution

Further information : Avoid subsoil penetration.

13. DISPOSAL CONSIDERATIONS

Product : Dispose of the product according to the defined EWC (European Waste Code) No.

Contaminated packaging : Take empty packaging to the recycling plant.

Waste key for the unused product : EWC 070601

Waste key for the unused product(Group) : Waste material of HZVA from fats, lubricants, soaps, detergents, disinfectants and personal protection products.

14. TRANSPORT INFORMATION

Further information : Not classified as dangerous in the meaning of transport regulations.

15. REGULATORY INFORMATION

General advice : In the EU, this product falls under the Directive medical devices 93/42/EEC. The product is classified and labelled in accordance with EC directives or respective national laws.

Symbol :



Xi
Irritant

R-phrases) : R36/38 Irritating to eyes and skin.

S-phrases) : S26 In case of contact with eyes, rinse immediately with plenty of water and seek medical advice.

S37 Wear suitable gloves.

Special labelling of certain preparations : Labelling according to Regulation (EC) No. 648/2004:., < 5 % non-ionic surfactants, < 5 % soap, perfumes

National legislation

Legislation on the control of major-accident hazards involving dangerous substances : Directive 96/82/EC does not apply

SAFETY DATA SHEET

according to Regulation (EC) No. 1907/2006

Aspirmatic **No Change Service!**

Version 01.07

Revision Date 23.03.2009

Print Date 02.09.2009

Volatile organic compounds : 1 %
(VOC) content Directive 1999/13/EC on the limitation of emissions of volatile organic compounds

VOC-Content : none, Switzerland. VOC Ordinance, Annex II (Products)

16. OTHER INFORMATION

Application : Information on uses can be found in a separate preparation / product brochure.

Text of R-phrases mentioned in Section 3

R11 Highly flammable.
R22 Harmful if swallowed.
R34 Causes burns.
R50 Very toxic to aquatic organisms.

The information given is only of use for possible safety requirements and is based on our actual knowledge. This information does not guarantee any properties. It is the responsibility of the receiver / user of our products to follow existing laws and regulations.

|| Changes compared with the previous edition!!!