

SAFETY DATA SHEET

according to Regulation (EC) No. 1907/2006

schülke -+

Mikrozid Sensitive Wipes

No Change Service!

Version 01.01

Revision Date 15.01.2008

Print Date 25.02.2009

1. IDENTIFICATION OF THE SUBSTANCE/PREPARATION AND OF THE COMPANY/UNDERTAKING

Product information

Trade name : Mikrozyd Sensitive Wipes
Producer/Supplier : Schülke & Mayr UK Ltd.
Cygnet House
1, Jenkin road
S9 1AT South Yorkshire
United Kingdom
Telephone: +441142543500
Telefax: +441142543501
mail.uk@schulke.com
www.schulke.com

Contact person : Product / Application advice
+49 (0)40/ 521 00 666
pab@schuelke.com

Emergency telephone : +441142543500
+49 (0)40 / 52 100 -0

Application : Disinfectants, Biocidal product

2. HAZARDS IDENTIFICATION

Risk advice to man and the environment

Contact with the eyes can cause irritation.

3. COMPOSITION/INFORMATION ON INGREDIENTS OF THE ACTIVE SOLUTION

Chemical nature : (Preparation)
Aqueous containing solution on cellulose wipes

No dangerous ingredients according to Regulation (EC) No. 1907/2006

4. FIRST AID MEASURES

Eye contact : In the case of contact with eyes, rinse immediately with plenty of water and seek medical advice.
Skin contact : Wash with water and soap as a precaution.

5. FIRE-FIGHTING MEASURES

Special protective equipment for fire-fighters : No special precautions required.
Suitable extinguishing media : Water, Dry powder, Carbon dioxide (CO₂), Foam

Mikrozyd Sensitive Wipes**No Change Service!**

Version 01.01

Revision Date 15.01.2008

Print Date 25.02.2009

- Extinguishing media which shall not be used for safety reasons : No information available.
- Specific risk from the substance or the product itself, its combustion products or evolved gases : No special risks to be expected.

6. ACCIDENTAL RELEASE MEASURES

- Personal precautions : No special precautions required.
- Environmental precautions : No special environmental precautions required.
- Methods for cleaning up : Use mechanical handling equipment.

7. HANDLING AND STORAGE**Handling**

- Safe handling advice : not required under normal use
- Advice on protection against fire and explosion : No special protective measures against fire required.

Storage

- Requirements for storage areas and containers : Store at room temperature in the original container.
- Further information : Keep container tightly closed. Keep away from heat. Keep away from direct sunlight.
- Advice on common storage : Keep away from food and drink.

8. EXPOSURE CONTROLS/PERSONAL PROTECTION**Components with workplace control parameters**

none

Personal protective equipment

- Hand protection : Prolonged contact: Nitrile rubber gloves e.g. Camatril (>480 Min., Layer thickness: 0,40 mm) or butyl rubber gloves e.g. Butoject (>480 Min., Layer thickness: 0,70 mm) made by KCL or gloves from other manufacturers offering the same protection.
- Hygiene measures : Keep away from food and drink.
- Protective measures : Avoid contact with eyes.

9. PHYSICAL AND CHEMICAL PROPERTIES OF THE ACTIVE SOLUTION**Appearance**

- Form : liquid
- Colour : colourless
- Odour : characteristic

Mikrozid Sensitive Wipes

Version 01.01

No Change Service!

Revision Date 15.01.2008

Print Date 25.02.2009

Other data

solidification range	:	ca. 0 °C
initial boiling point	:	ca. 100 °C
Flash point	:	not applicable
Density	:	ca.0,998 g/cm ³ at 20 °C
Water solubility	:	at 20 °C in all proportions
pH	:	6 - 8 at 1.000 g/l (20 °C)

10. STABILITY AND REACTIVITY

Hazardous reactions	:	None reasonably foreseeable.
Hazardous decomposition products	:	None reasonably foreseeable.

11. TOXICOLOGICAL INFORMATION OF THE ACTIVE SOLUTION

Further information	:	No data is available on the product itself. The classification was made according to the calculation procedure of the Preparations Directive.
---------------------	---	---

12. ECOLOGICAL INFORMATION OF THE ACTIVE SOLUTION

Further information	:	This product has no known eco-toxicological effects. The surfactant(s) contained in this preparation complies(comply) with the biodegradability criteria as laid down in Regulation (EC) No.648/2004 on detergents. Data to support this assertion are held at the disposal of the competent authorities of the Member States and will be made available to them, at their direct request or at the request of a detergent manufacturer. Avoid subsoil penetration. Do not flush into surface water or sanitary sewer system.
---------------------	---	---

13. DISPOSAL CONSIDERATIONS

Product	:	Can be incinerated or landfilled together with household waste in compliance with the regulations, and after consultation with the waste disposal services.
Contaminated packaging	:	Take empty packaging to the recycling plant.

Mikrozid Sensitive Wipes

Version 01.01

No Change Service!

Revision Date 15.01.2008

Print Date 25.02.2009

14. TRANSPORT INFORMATION

Further information : Not classified as dangerous in the meaning of transport regulations.

15. REGULATORY INFORMATION

General advice : The product does not need to be labelled in accordance with EC directives or respective national laws.

National legislation

Legislation on the control of major-accident hazards involving dangerous substances : Directive 96/82/EC does not apply

VOC-Content : none, Directive 1999/13 /EC on the limitation of emissions of volatile organic compounds

VOC-Content : none, Switzerland. VOC Ordinance, Annex II (Products)

16. OTHER INFORMATION

The information given is only of use for possible safety requirements and is based on our actual knowledge. This information does not guarantee any properties. It is the responsibility of the receiver / user of our products to follow existing laws and regulations.

|| Changes compared with the previous edition!!!